

Solid Phase Synthesis of Peptoids

GenScript USA Inc.

Peptide drugs are a rapidly growing class of therapeutics. However, the peptide drug discovery was hampered by its inherent characteristics: low stability to enzymatic digestion, high conformational flexibility resulting in low specificity to target, low hydrophobicity and lack of specific transportation systems.

In order to overcome the above disadvantages of peptides, peptidomimetics are driving extensive research efforts and have shown a promising prospect in peptide drug discovery. A peptidomimetic is a compound designed to mimic the activity of a peptide, such as a peptoid and a beta-peptide. But a peptidomimetic should have structural difference from its parent peptide. This structural difference can bring us greater advantages for its function as a drug. Many classes of peptidomimetics have been reported [for review, see ref. 1, 2].

Among these peptidomimetics, peptoids are completely resistant to proteolysis and easily get a great number of diversity. Peptoids are oligomers N-substituted glycines developed by Zuckerman and colleagues in 1992 [ref.3]. Peptoids represent a new class of synthetic compound. Comparing with a peptide, as we can see in Figure 1, the side chains are connected to the nitrogen of a peptide backbone in a peptoid instead of alpha carbon in a peptide. In native peptides, R represents 20 kinds of substitution, while in peptoids, R is roughly unlimited.




Figure 1: Structural comparison between a peptide and a peptoid

The synthesis of a peptoid can be achieved through monomer and submonomer methods. The monomer method is exactly the same as solid phase peptide synthesis except that all the protected monomers for peptoid synthesis need to be synthesized first. Zuckerman and colleagues provided three routes illustrated in Figure.2 [ref.3] to prepare peptoid monomers. The synthesis of a peptoid oligomer was performed using the standard Fmoc peptide synthesis method with in situ activation by PyBOP or PyBroP. Since all the backbone nitrogens are alkylated, the coupling is more difficult than normal peptide synthesis. That is why PyBOP instead of conventional carbodiimide coupling reagents was used. Using this method, ten peptoid oligomers were synthesized and showed similar activity with their parent peptides and higher protease resistance [ref.3]. The monomer method is quite straightforward. However, the synthesis of every N^α-protected monomer is tedious and time-consuming.




Figure 2: General route for synthesis of peptoid monomers

The submonomer method [ref.4] eliminates the need of N^z-protected monomers. It consists of two step simple reactions: an acylation step performed by the addition of bromoacetic acid and DIC and a nucleophilic displacement of bromide with a primary amine as can be seen in Figure 3. Bromoacetic acid and primary amines are readily available submonomers. Because this approach shows greater advantages than the monomer approach, it has been widely applied to peptoid library construction [ref. 5, 6].





Figure 3: Submonomer route to synthesize peptoids

Both of the above methods are quite straightforward and easily adapted to automated synthesis. However, the diketopiperazine formation and reaction of stereo-hindered monomers or submonomers are still problematic. Microwave technology was applied to accelerate the reaction especially for the hindered amines [ref. 7, 8]. Selection of a stereo-hindered resin such as a 2-chloro-trityl resin can effectively avoid the formation of diketopiperazine during peptoid synthesis. Microwave technology combined with a multi-channel peptide synthesizer enables us get hundreds of peptoids within hours. The “Split-mix” method (for review, see ref. 9) further generates large combinatorial libraries of peptoids.

As mentioned above, peptoids show enhanced proteolytic stabilities and increased cellular permeabilities relative to their parent peptides. With the progress of peptoid high-throughput synthesis and screening, peptoids hold significant promises for therapeutic applications (for review, see ref. 10).

References

- (1) Vagner, J. et al, *Curr. Opin. Chem. Biol.* **2008**, 12(3), 292.
- (2) Soth, M. J. and Nowick, J. S. *Curr. Opin. Chem. Biol.* **1997**, 1, 120.
- (3) Simon R. J.; et al, *Proc. Natl. Acad. Sci. U.S.A.* **1992**, 89, 9367.
- (4) Zuckermann, R. N. et al, *Journal of the American Chemical Society*, 114(26), 10646.
- (5) Zuckermann, R. N., et al, *J. Med. Chem.* **1994**, 37, 2678.
- (6) Martin, E.J., et al, *J. Med. Chem.* **1995**, 38, 1431.
- (7) Hernando J. et al, *Org. Lett.*, **2002**, 4 (23), 4057.
- (8) Benjamin C. et al, *Org. Lett.*, **2005**, 7 (8), 1521
- (9) Lam, K. et al, *Chem. Rev.*, **1997**, 97 (2), 411
- (10) Sarah, A. et al, *Org. Biomol. Chem.*, **2009**, 7, 1508-1524



Toll-Free: 1-877-436-7274

Tel: 1-732-885-9188

Fax: 1-732-210-0262

Email: peptide@genscript.com

Web: www.genscript.com

GenScript USA Inc. is a leading biology CRO focusing exclusively on early drug discovery and development services. Built on our assembly-line mode, one-stop solution, continuous improvement, and stringent IP protection, GenScript provides a comprehensive portfolio of services that include Bio-Reagent, Bio-Assay, Lead Optimization, and Antibody Drug Development. GenScript has over 1,000 employees and is headquartered in New Jersey, USA with a major subsidiary in China. For more information about GenScript USA Inc., visit: www.genscript.com

Custom Peptide Services

—FlexPeptide™ Ensures Delivery & Quality!

GenScript's proprietary FlexPeptide™ platform provides custom peptides of up to 200 residues in length and a capacity of 6,000 peptides per month. With a flexible scale ranging from milligrams to kilograms, comprehensive labeling, modification options, and rapid high-throughput peptide library synthesis, this platform is your ideal choice for synthesizing peptides and/or a peptide library. GenScript also offers large-scale cosmetic peptide and cGMP-grade production services.

We, at GenScript, provide high quality peptide synthesis service with a success rate well above the industry standard because we are committed to Total Quality Management (TQM) to assure customers' complete satisfaction.

Peptide Synthesis

—Fully quality assured with guaranteed delivery!

- Long Peptides Up to 200 AA
- Instant online quote and order
- Free peptide antigen design

Peptide Modification

—Comprehensive modifications with competitive prices

- Amidation and Acetylation
- Biotin and FITC labeling
- PEGylation
- KLH, BSA, OVA conjugation, etc.

Peptide Library Services

—One-Stop solutions for your peptide-based drug research!

- Flexible purity and quantity choices
- Instant online quote
- One-Stop down-stream services including epitope mapping and assay

Peptide Array Services

—One-Stop solutions for mapping protein interactions!

- SPOT technology
- No cross-contamination
- Stringent quality control

cGMP Peptide Synthesis

—Individual solutions for projects of any size!

- Advanced FlexPeptide™ technology
- Proven track records
- Customized solutions

Click Peptide Service


—Turn hydrophobic to hydrophilic by one “CLICK”!

- Increase the solubility of hydrophobic peptide
- Slow the aggregation of β -amyloid peptide
- Applicable to cell signal transduction and other studies





GenScript's Total Quality Management (TQM) System for Peptide Services



Long peptide case study

The FlexPeptide™ system employs GenScript's proprietary ligation method to achieve high yields even for extremely long peptides. This is accomplished by, first synthesizing several shorter sequences and then efficiently ligating them together.


Sequence: Ala-Glu...-Lys-...-Lys-Lys-...-Glu-Ala-...-Ala-Ala

Length: 149 AA


M.W.: 15,971.51

Purity: 93%

HPLC results:



MS results:



For more information on Peptide Services, visit http://www.genscript.com/peptide_overview.html

